

Classic Scripting in Pinnacle³

Michael Meltsner, Karl Bzdusek, Matthieu Bal, Ying Xiong

Philips Radiation Oncology Systems

May 7, 2014

Version 1.0

PHILIPS

Disclaimer

For non-clinical use

- The following presentation outlines advanced scripting techniques
- Philips officially only supports “Record and Playback” within clinical releases of Pinnacle
- The tools that follow are meant for research mode only and are not validated for use on patients

Overview

What is a Script?

- Scripting in Pinnacle³ gives the user the ability to communicate to/from the Pinnacle core without having access to the code directly
- Each script is an ASCII file that contains a collection of messages that communicate sequentially with Pinnacle³
- These messages allow the user to mimic interactions with mouse clicks and keyboard input
- They are analogous to macro type functions in other software like Excel

Scripting Messages


Three basic types

- Pinnacle responds to 3 basic types of Scripting Messages:
 - 1) **Set command:** tells Pinnacle³ to change an underlying attribute
 - Beam angle, prescription name
 - 2) **Query request:** asks Pinnacle³ to return the value of an attribute
 - Number of beams, Patient's name
 - 3) **Action command:** tells Pinnacle³ to perform a task
 - Add a beam, Delete an ROI

Pinnacle Architecture

Pinnacle³ MVC Architecture


Set, Query, and Action messages



Pinnacle³ GUI

How widgets work

- The Pinnacle³ GUI is essentially an extensive set of messages hidden behind widgets
- Each widget in the GUI has 1 or more of the three types of messages associated with it
- The label widget displays the result of a **Query** request
- The button widget in Pinnacle³ typically performs an **Action**
- The text box displays a **Query** request as its value and sends a **Set** message after the user enters a new value



Pinnacle³ Internal Design

Understanding scripting messages




- Pinnacle³ is designed in an object-oriented manner
- Every object has associated attributes and functions
- The hierarchy of the Pinnacle³ objects is important and directly relates to the functionality of the software and hence is mimicked in the scripting language
- Objects of the same class are stored in a container simply called **ObjectList**.
 - Each object within that list may contain subobjects including sublists

Pinnacle³ Internal Design

Hierarchy

- The hierarchy of Pinnacle³ (and thus Scripting) follows a logical sense of design
 - For example, each plan may contain multiple **Trials**. Each **Trial** may contain multiple **Beams**. Each **Beam** may contain multiple **MLC Leaves** and so on.



- You can **Set, Query, or Act** on any object by following the hierarchical structure of the Pinnacle³ design starting from the root level, in this case **TrialList**.

Pinnacle³ Internal Design

Root level objects

- The top most objects are called “Root level”
 - Examples of Root Level **ObjectLists**
 - TrialList - Root level container of trials
 - RoiList - Root level container of ROIs
 - DVHList - Root level container of DVH
 - VolumeList - Root level container of volumes
 - ViewWindowList - Root level container of view windows
 - WindowList - Root level container of windows
- An example of a non-root container follows from the previous slide
 - TrialList.IMRT.BeamList ...- Each **Trial** contains a **BeamList**
- **Note: Not all Root level objects are ObjectLists**
 - PlanInfo - Root level object that holds information about the plan/patient

Message Syntax and Examples

Basic Syntax

Scripting rules

- Each message within the script must end with a semi colon
- With some exceptions, there is one “=” per line for each message
- All script commands are case sensitive
- The script may contain comments denoted by starting a line with a “//”
`//This is my first script!`

Basic Syntax

Set message

- The set message sends information to Pinnacle³
- It contains the attribute you want to set on the left and its new value on the right

<SetMessage> = <Value>;

<Value> can be one of the following:

- A string enclosed in double quotes
 - A floating point or integer number
 - A Query message
-
- **Set** the color of the ROI named “Bladder” to green
RoiList.Bladder.Color = “green”;

Basic Syntax

Query request

- The query request asks Pinnacle³ to return some information about an object
- The query request is most commonly used as input for a **Set** message

<SetMessage> = <Query Request>;

- **Set** gantry of “Beam_1” equal to gantry of beam “AP” via **Query**
TrialList.IMRT.BeamList.Beam_1.Gantry =
TrialList.IMRT.BeamList.AP.Gantry;

Basic Syntax

Action message

- The **Action** command tells Pinnacle³ to perform a task
- The **Action** message ignores the right side of the equals sign, but still requires one to properly execute

<ActionMessage> = "";

where *"" is just an empty string indicated by quotations*

- Recompute the statistics of an ROI via **Action message**
RoiList.Current.RecomputeStatistics = "";

Example Messages

Reference specific objects

- Special name **Current**, **First**, **Last** messages highlighted object in GUI's list

TrialList.Current.BeamList.Last.Couch = 180.0;

- Reference by name (delimiter necessary if name has a space)

TrialList.Current.BeamList.Lateral.Couch = 180.0;

OR

TrialList.Current.BeamList.# "Lateral".Couch = 180.0;

TrialList.Current.BeamList.# "Ant Post".Couch = 180.0;

- Reference the 3rd element in the list by Index

TrialList.Current.BeamList.# "#2".Couch = 180.0;

- Special character * loops over all objects in list

TrialList.Current.BeamList.# "*".Couch = 180.0;

TrialList.# "*".BeamList.# "*".Couch = 180.0;

Creating Children

Expand the Objectlist

- A **child** is an object in a list (e.g. **Prescription** in a **PrescriptionList**)
- The **CreateChild Action** message appends a new instance of the child class to the list
TrialList.Current.PrescriptionList.CreateChild = "";
TrialList.Current.PrescriptionList.CreateChildMakeCurrent = "";
- Alternatively, send the class name to the **ObjectList**
TrialList.Current.PrescriptionList.Prescription= "";
- Exception: the **CreateChild Actions** for adding beams and regions of interest
CreateNewBeam = ""; **CreateNewROI = "";**

Remove Children

Contract ObjectList

- Destroy the current child and remove it from the list
`TrialList.Current.BeamList.Current.Destroy = "";`
- Remove the current child from the list (don't free memory)
`TrialList.Current.BeamList.Current.Remove = "";`
- Make the last child current
`TrialList.Current.BeamList.Last.MakeCurrent = "";`
- Print to terminal the number of elements in the list
`Echo = TrialList.Current.BeamList.Count;`
- Add 3 children (Note: this won't work for beams or ROI's)
`PoiList.AddChildren = 3;`

Sorting Children

Order the objects intelligently



- An **ObjectList** can be sorted by any child object key
- Sort the beam list by couch angle in ascending order
TrialList.Current.BeamList.SortBy.Couch = "";
- Multiple keys can be specified
TrialList.Current.BeamList.SortBy.Couch.Gantry = "";
- To do a descending sort, precede the attribute with "D"
TrialList.Current.BeamList.SortBy.Couch.D.Gantry = "";
- The primary sort key is the couch angle in ascending order
- The secondary sort key is the gantry angle in descending order
- All capitals in names are higher in sort order than lowercase

Root Level Commands

Examples

- Pinnacle³ contains various root level commands
- These commands tell Pinnacle³ to perform some **Action** on the entire plan

- Exiting the System

```
Quit = "";  
QuitWithSave = "";
```

- Saving Plan

```
SavePlan = "";
```

- Echo String to the Console (i.e. print statement)

```
Echo = "Starting My Script";
```

GUI Control Messages

Examples

- Displaying Wait Messages in Status Bar
WaitMessage = "Doing something slow...";
...
WaitMessage = "Doing something slower...";
...
WaitMessageOff = "";
- Issuing a Warning Message (Requires user interaction)
WarningMessage = "This is a test.";
- Issuing an AskYesNo Window
AskYesNoPrompt = "Are you having fun?";
AskYesNoDefault = 1; // Default yes
IF.AskYesNo.THEN.WarningMessage = "Having fun!";

Dismissing Warnings and Yes/No messages

Examples

- Dismiss the next yes/no message and reply yes
Test.ExpectAskYesNo = 1;
Test.ExpectedAskYesNoReply = 1; // Yes
- Dismiss the next yes/no message with the text “Warning text” and reply no
Test.ExpectAskYesNo = “Warning text”;
Test.ExpectedAskYesNoReply = 0; // No
- Dismiss the next warning message
Test.ExpectWarningMessage= 1;
- Dismiss the next warning message with the text “Warning text”
Test.ExpectWarningMessage= “Warning text”;

Creating Your Own Variables

Using the “Store”


Creating your own variables

- Pinnacle³ allows the user to create 3 types of variables that can be used inside your script. These three types are:
 - 1) **Float** - numerical or boolean operations
 - 2) **String** - character based information to be stored and manipulated
 - 3) **Pointer** - references something in Pinnacle³ memory directly
- These three variables are collectively maintained in the Pinnacle³ “**Store**”
- The **Store** is not saved with the plan by default
- The **Store** can be a root level Object, or can be created for specific, existing Pinnacle³ Objects

Store Variables

Float creation/deletion

- Creating Float Entry
Store.FloatAt.MyFloat = 1.23;
- Modifying Float Entry
Store.At.MyFloat.Value = 2.34;
OR **Store.FloatAt.MyFloat = 2.34;**
- Deleting a Float Entry in the store
Store.FreeAt.MyFloat = "";
- Printing a float value to the xterm
Echo = Store.At.MyFloat.Value;



Store Variables

Float numeric operations

- Add, Subtract, Multiply, and Divide syntax
Store.At.MyFloat.Add = 23.0; //Adds 23 to MyFloat
- Square, SquareRoot, Negate, and Invert(1/x) Syntax
Store.At.MyFloat.SquareRoot = ""; //Replaces MyFloat with SQRT(MyFloat)
- Round, Nint, and Absolute
Store.At.MyFloat.Round = ""; //Rounds MyFloat
- Can also be used as a Query
Store.FloatAt.MyFloat = Store. At.MyOtherFloat.Round;
- Query attributes of Pinnacle³ objects to set the value of the float
Store.FloatAt.MyFloat = TrialList.Current.BeamList.Current.Gantry;

Store Variables

String operations

- Creating a String Entry

Store.StringAt.MyString = "ABC";

- Modifying a String Entry

Store.At.MyString.String = "DEF";

OR Store.StringAt.MyString = "DEF";

- Concatenating Strings

Store.At.MyString.AppendString = RoiList.Current.Name;

- Deleting a String Entry

Store.FreeAt.MyString = "";

- Getting Environment Variables

Store.StringAt.HomeDir = GetEnv.HOME;

Store Variables

Pointers

- Creating the Reference by querying the memory address of an Object

```
Store.At.MyBeamReference =  
    TrialList.Current.BeamList.Current.Address;
```

- Modifying the Reference

```
Store.At.MyBeamReference.Gantry = 123;  
Store.At.MyBeamReference.Couch = 234;
```

- Removing the Reference

```
Store.RemoveAt.MyBeamReference = "";
```

- **Note: Do not use FreeAt; it destroys the Beam object**

Advanced Techniques

Conditional Operations

Structure

- **IF.key.THEN.action.ELSE.message = value;**
Where key is a boolean query.
- **IF.key.THEN.message = value;**
- **IF.!key.THEN.message = value;**
Where ! is used for to check for false boolean keys



Comparator Operations

Compare values

- **IF.key1.OP.key2.THEN.action.ELSE.message = value;**
- “OP” can be:

Numeric	String	Logical
EQUALTO	Is	AND
NOTEQUALTO	IsNull	OR
GREATERTHAN	Contains	XOR
LESSTHAN	STRINGEQUALTO	
GREATERTHANOREQUALTO	STRINGNOTEQUALTO	
LESSTHANOREQUALTO		

Conditional and Comparator

Examples

- Is:

IF.Store.At.MyString.Is.# "xyz".THEN.WarningMessage = "It's xyz!";

- Contains:

IF.Store.At.MyString.Contains.# "xyz".THEN.WarningMessage = "Has xyz!";

- **Note: Right side must be a string not a variable**
- Use **STRINGEQUALTO** for two queries instead of Is

- IsNull:

IF.Store.At.MyString.IsNull.THEN.WarningMessage = "It's empty!";

- IF - THEN - ELSE numeric:

**IF.TrialList.Count.EQUALTO.# "#1".THEN.
TrialList.CreateChild.ELSE.TrialList.Last.MakeCurrent=1;**

Executing Scripts from Within Scripts

Examples

- Run another script via message

```
ExecuteNow = "/home/pinnbeta/MyScript.Script";
```

- Note the use of two “=” signs

```
Store.StringAt.MyCommand =  
    "RoiList.Current.Density = RoiList.Bone.Density";  
Store.At.MyCommand.Execute = "";
```

- Use “#” to delimit a string within a string

```
Store.StringAt.MyCommand = "RoiList.Current.Name = #Bone";  
Store.At.MyCommand.Execute = "";
```

- Use **LoadNoChecksum** with a string variable if you want to execute another script

```
Store.StringAt.MyCommand =  
    "LoadNoChecksum = /usr/tmp/AddToCount.Script";  
Store.At.MyCommand.Execute = "";
```


Edit the Store Via the GUI

Ask user for input Via window

- Create a window to display and edit variables

WindowList.MyStoreEditor.CreateStoreEditor = "Store.Address";

- Supports Strings and Floats
- Re-launch using the same command
(existing window will be created)
- Parameters can be edited by the user.
- Script continues after window launches
without waiting for user input



Force User Entry into the Store Editor

Stop Script from continuing until entry completed

- Create a window to display and edit variables

WindowList.MyStoreEditor.CreateStoreEditor = "Store.Address";

- Dismiss the window

WindowList.MyStoreEditor.Unrealize = "";

- Make it modal

WindowList.MyStoreEditor.IsModal = 1;

- Re-launch the window

WindowList.MyStoreEditor.Create = "";

...

// User dismisses the store window and script continues

Add Attributes to Existing Objects

Append custom variables to any Pinnacle object

- The custom **Store** is specific to each object, not to the class (e.g. two beams might have different variables)
- Functionality is identical to the root **Store**, except that the data is persisted (i.e. saved and copied with object)
- Create an attribute specific to the current beam called MyFloat with value 123.4
**TrialList.Current.BeamList.Current.
Store.FloatAt.MyFloatAttribute = 123.4;**
- You may add as many custom attributes as desired to an existing Pinnacle³ object and edit them via the **Store Editor**

Adding a Dependency

Trigger something to happen automatically

- Execute MyRoiScript.Script whenever the number of ROIs changes

```
KeyDependencyList.CreateChild = "";
```

```
KeyDependencyList.Last .Name = "NewRoiDependency";
```

```
KeyDependencyList.Last .KeyString = "RoiList.Count";
```

```
KeyDependencyList.Last .AddAction =
```

```
    "Script.ExecuteNow = #/home/pinnbeta/MyRoiScript.Script";
```

Iterating

Method 1: using the * and @ characters



- Special character * loops over all objects in list setting one attribute = 180.0

TrialList.Current.BeamList.# "*" .Couch = 180.0;

- Make each ROI current, then execute another script using @

**RoiList.ChildrenEachCurrent.# "@" .Script.ExecuteNow =
"/home/pinnbeta/MyScript.Script";**

Iterating

Method 2: use recursive scripts part 1

- **RecursiveLoopExample.Script** calls iteratively the **AddToCount.Script**
**Store.StringAt.ScriptCommand = "LoadNoCheckSum =
/home/pinnbeta/source/Scripts/AddToCount.Script";**
- Initialize the counter and limit
Store.FloatAt.Counter = 0.0;
Store.FloatAt.Limit = 10.0;
- Start the loop
Store.At.ScriptCommand.Execute = "";
- Free the store variables
Store.FreeAt.ScriptCommand = "";
Store.FreeAt.Counter = "";
Store.FreeAt.Limit = "";

Iterating

Method 2: Use recursive scripts part 2

- *AddToCount.Script*

```
Store.At.Counter.Add = 1.0;
```

```
Echo = "Adding to counter...";
```

```
Echo = Store.FloatAt.Counter;
```

```
IF.Store.FloatAt.Counter.LESSTHAN.Store.FloatAt.Limit.
```

```
THEN.Store.At.ScriptCommand.Execute = "";
```


Initialization Scripts

Have scripts run every time Pinnacle³ opens

- There are initialization scripts that can be created to run at Pinnacle³ and Launchpad startup.
 - **LaunchpadInit** - all commands are run after the database is loaded
 - **PinnacleInit** - all commands in the script are run after the plan is loaded
- These scripts do not exist by default
- Create them and put them in your home directory e.g. /home/pinnbeta

Communicating with the Operating System

Examples

- Write disk usage info to a text file
SpawnCommand = "df -k > /usr/tmp/MyTempFile.txt";
- Write text to a text file (overwrites contents) of file
SpawnCommand = "echo HelloWorld > /usr/tmp/MyTempFile.txt";
- Write text to a text file (appends contents)
SpawnCommand = "echo HiAgain >> /usr/tmp/MyTempFile.txt";
- Launch the gedit text editor and view file as background
SpawnCommand = "gedit /usr/tmp/MyTempFile.txt &";
- Execute command in terminal while Pinnacle continues to run
SpawnCommandNoWait = "xterm -e OSCommand";

Writing Out Specific Data to OS

Write select Pinnacle³ data to text file via OS


- Find the DVH's first data point and write it out to text file

```
Store.StringAt.EchoCommand = "echo ";  
Store.FloatAt.Dose = DVHList.Current.Data.PointsFirst;  
Store.At.EchoCommand.AppendString = Store.At.Dose.Value  
Store.At.EchoCommand.AppendString = " >>  
    /usr/tmp/OutputFile.txt";  
SpawnCommand = Store.At.EchoCommand.String;
```

How to Run and Create Your Own Scripts

How to Run Your Script


Opening the scripting window



Finding Messages

Record and playback your GUI interactions

- Recording a script writes out to a text file all the messages that Pinnacle³ is executing via user interaction
- You may view, edit, rerun the recorded script
- Recording scripts is an easy way to learn the Pinnacle³ scripting language and find specific messages



PHILIPS

Finding Messages

Read/Write objects to ASCII files on disk

- Writing out an object will place all the persisted data into a text file
`TrialList.Current.BeamList.Current.Save = "/files/rtp/Beam.out";`
- Read all the data back in from the text file
`TrialList.Current.BeamList.Current.Load = "/files/rtp/Beam.out";`
- Save and Load the custom **Store**
`Store.Save = "/home/pinnbeta/MyStore.Store";`
`Store.Load = "/home/pinnbeta/MyStore.Store";`
- **Note: The existing Store is not emptied before loading and variables are overwritten if identical names exist**

Finding Messages

Examine Pinnacle³ files




- Nearly all data within a Pinnacle³ plan are stored to disk
- Pinnacle³ stores nearly all of its data in Script format, i.e. a series of scripting messages exactly like the ones you would create on your own
- You may examine any of these files to find all the persisted objects, attributes, and functions that can be messaged
- These files include:
 - **Transcript Files** - /Institution_X/Mount_Y/Patient_Z/Plan
 - **Plan Files** - Institution_x/Mount_y/Patient_z/Plan_a/plan.Trial
 - **GUI Configuration files** - PinnacleStatic/Config/*.cfg
- **Note: Do not modify these files!**

Finding Messages

Using the ToolTips

- Pinnacle³ can display to you the scripting messages behind every widget in the GUI simply by hovering your mouse
- Enable **ToolTips** in Planning mode (available only in research version after 3/2014)
ShowToolTips = 1;
- To enable this in Launchpad, add the same message to the **LaunchpadInit** file



Finding Messages

Search reference materials

- As part of this training, your Philips instructor provides additional documentation such as:
 - **Advanced scripting manual**
 - **Example scripts**
 - **Text files with lists of the most common objects and all their messages**



THANK YOU